

# Ventilators

# ZOLL®



## ACCESSORIES AND DISPOSABLES

# ACCESSORIES FOR ALL MARKETS

Ensures safe and effective treatments for infants through adults



## CPAP MASK AND HARNESS

Interface between the circuit and the patient for noninvasive ventilation on CPAP or BL (CPAP + PS) modes. The mask is nonvented. It includes the harness.

CPAP Harness + Mask #5  
Regular Adult  
712-0002-00  
712-0002-20 (case of 20)  
712-0002-50 (case of 50)

CPAP Harness + Mask #6  
Large Adult  
712-0003-00  
712-0003-20 (case of 20)  
712-0003-50 (case of 50)

CPAP Harness + Mask #4  
Child  
712-0004-00  
712-0004-20 (case of 20)  
712-0004-50 (case of 50)

## CPAP MASK

Interface between the circuit and the patient for noninvasive ventilation on CPAP or BL (CPAP + PS) modes. The mask is nonvented.

CPAP Mask #1 Small Infant  
812-0006-00  
812-0006-20 (case of 20)


CPAP Mask #2 Infant  
812-0007-00  
812-0007-20 (case of 20)

CPAP Mask #3 Small Child  
812-0008-00  
812-0008-20 (case of 20)

CPAP Mask #4 Child  
812-0009-00  
812-0009-20 (case of 20)

CPAP Mask #5 Regular Adult  
812-0010-00  
812-0010-20 (case of 20)

CPAP Mask #6 Large Adult  
812-0011-00  
812-0011-20 (case of 20)



## CPAP HARNESS

The CPAP harness is the strap used to attach the CPAP mask to the patient for hands-free noninvasive ventilation.

CPAP Harness  
334-0125-00  
334-0125-10 (case of 10)

## BV FILTER

50 bacterial viral filters. May be used on the Fresh Gas intake.

Bacterial/Viral Filter  
465-0024-00 (case of 50)



## AIR INLET FILTER

Internal disk filter, needs to be changed when a Fresh Gas intake alarm is present.

Internal Disk Filter  
465-0027-00  
465-0027-50 (case of 50)



## VENTILATOR CIRCUITS

Six-foot disposable, single patient use ventilator circuits.

Adult Circuits  
820-0106-00  
820-0106-15 (case of 15)



Infant/Pediatric Circuits  
820-0107-00  
820-0107-20 (case of 20)



## HEAT MOISTURE EXCHANGERS

Heat moist exchanger with filter; provides humidification to the patient during ventilation. Not for use with CPAP or BL.

HME Adult  
465-0025-00 (case of 20)  
820-0053-00 (case of 50)



HME Pediatric Dead Space <25ML  
820-0109-00 (case of 20)



HME Adult Dead Space <75ML  
820-0108-00 (case of 20)



HME Infant Dead Space <10ML  
820-0110-00 (case of 30)



## PULSE OXIMETER CABLES

Masimo Pulse oximeter sensors to use with the ventilator.

LNCS Adult Reusable  
SpO<sub>2</sub> Sensor  
8000-0294



LNCS Disposable Adult SpO<sub>2</sub>  
Sensors (Case of 20)  
8000-0320



LNCS Pediatric Reusable  
SpO<sub>2</sub> Sensor  
8000-0295



LNCS Pediatric Disposable  
SpO<sub>2</sub> Sensors (Case of 20)  
8000-0321



LNCS Inf-3, Infant Disposable SpO<sub>2</sub>  
Adhesive Sensor (Case of 20)  
8000-0322



Four-foot LNCS Reusable  
SpO<sub>2</sub> Patient Cable  
8000-0298



## OTHER ACCESSORIES

Other standard accessories for ZOLL ventilators

### POWER SUPPLY

Power supply and power cords for ZOLL ventilators.

Power Supply  
024-0012-00



Power Cord  
708-0042-00



Power Supply  
and Power Cord  
710-0731-01



### OTHER

Test lung for demonstrations and product education.

One-Liter Dark Blue Test Lung  
820-0132-00



Six-foot Oxygen Hose for  
High Pressure FiO<sub>2</sub> delivery  
825-0002-00



### CASES

Hard cases for ventilators and storage.

Ventilator Hard Case  
703-0731-03




Ventilator Hard Case with Wheels  
703-0731-02



# ACCESSORIES FOR HOSPITALS

Supports safe operation within the hospital



## CARRIER

White cage ventilator carrier for Eagle II.  
703-0731-27

## MRI CIRCUITS

12-foot Disposable MRI length circuits for ZOLL ventilators.

Adult Circuits  
820-0130-10 (case of 10)



Infant/Pediatric Circuits  
820-0131-10 (case of 10)



## VENTILATOR STANDS AND ACCESSORIES

Rolling cart used in hospitals as an accessory for MRI and non-MRI ventilators.

Circuit Holder  
820-0124-00  
(Used with non-MRI rolling cart 816-0731-00)



IV Pole  
704-0731-09  
(Used with rolling carts 816-0731-01 and 816-0731-00)



Rolling Cart  
816-0731-00



Rolling Cart (for MRI use)  
816-0731-01



# ACCESSORIES FOR EMS

Essential accessories for the EMS environment

## MOUNTING

Universal wall mount for ambulances  
(for US customers)  
703-0003-00



## CARRIER

Olive carrier cage for ventilators  
703-0731-17




## POWER SUPPLY

Vehicle charger  
704-0EMV-06



## LOW FLOW OXYGEN SYSTEM

Three-liter O<sub>2</sub> reservoir kit  
704-0004-00



# WHY THE ZOLL WEBSTORE?

Our WebStore offers convenient, one-stop shopping. You can view, check pricing, and order accessories and disposables for all of our product lines.

## FEATURES

### MY DEVICES

"My Devices" ensures you are ordering accessories that are compatible with the devices you own. Once set up, items will be flagged to indicate that they are compatible with your device or you will receive a warning message at checkout if not compatible.

### CONTRACT PRICING

Once logged in, your contracted price will be the price displayed.

### CREATE A QUOTE

Simply build a cart of items you want included in your quote. During the checkout process, a "Print Quote" button will appear that will generate a PDF of this information that you can print or email.

### CHECK ORDER STATUS

Easily view your order status, track shipping, and reorder with the click of a button.

### FREE SHIPPING

On WebStore orders over \$800



**Start Today!**  
[www.zollwebstore.com](http://www.zollwebstore.com)

## ZOLL MEDICAL CORPORATION

An Asahi Kasei Group Company | 269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | [www.zoll.com](http://www.zoll.com)

Copyright © 2017 ZOLL Medical Corporation. All rights reserved. AEV, Eagle II, EMV+ and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.  
PN 335

For subsidiary addresses and fax numbers, as well as other global locations, please go to [www.zoll.com/contacts](http://www.zoll.com/contacts).

# ZOLL®